



SOP – Tissue for cells

In case of an opportunity, from any chILD patient carrying a defined molecular mutation or not, fresh lung tissue or other tissue should be collected for the preparation of epithelial cells, and if not possible for fibroblasts for long term storage.

Opportunities: Operations, e.g. tracheostomy, PEG placement, etc
=> Skin / fibroblasts

Operations, e.g. Ductus ligation, lung biopsy, etc
=> lung / epithelial cells

Autopsy, selective lung only autopsy
=> lung / epithelial cells

Requirement: Written informed consent, e.g. chILD-EU form signed

Procedure: Put fresh, sterile tissue into culture medium (MEM; DMEM, RPMI, etc with antibiotics if available) and ship as fast as possible to Munich Biobank.

Get into contact **to arrange further processing of the samples**, with

Ms. Schams: +49 89 4400 53715, Andrea.Schams@med.uni-muenchen.de

Ms. Wesselak: +49 89 4400 57887, Traudl.Wesselak@med.uni-muenchen.de

Prof Griese: Matthias.griese@med.uni-muenchen.de

Further information and help, and sample shipping sheet to use: www.childeu.net